



LABORATORY KNOWLEDGE. DELIVERED.™

LABVANTAGE® REPORTING & ANALYTICS

LABVANTAGE® REPORTING & ANALYTICS

LABVANTAGE REPORTING & ANALYTICS

Today, enterprises are striving to significantly accelerate the decision-making process and reduce time-to-market by leveraging robust reporting and analytics tools. In information-rich industries, laboratories are under increasing pressure to deliver meaningful data with critical scientific and business insights. Real-time reporting and analytics provide scientific researchers, analysts, informaticians, project and business-line managers, and executive decision-makers the means to visualize data and more readily identify the critical issues and trends, ultimately enabling them to turn information into knowledge and make earlier, more informed critical decisions.

However, many traditional laboratory information management systems (LIMS) supporting laboratory operations do not offer direct, easy access to all the reporting and analytics functionality needed by both researchers and analysts. Lack of adequate access tends to hinder rather than promote information sharing and decision making in the organization. Without readily accessible reporting and analytics options, a users' ability to flexibly retrieve, distribute and manage valuable knowledge derived from laboratory data is restricted. Users have to deal with a multitude of disparate applications and analyze data with considerable manual effort, significantly impeding their ability to make real-time, informed decisions.

To optimize data accessibility and improve knowledge sharing, laboratories need a variety of integrated business intelligence tools. Users should be able to easily access, compile and analyze information via any number of reporting and analytics tools, as preferred, depending on the actual needs of the various management, operation, and research users and their requirements. With direct and easily accessible functionality, users can rapidly retrieve and combine data, enabling them to develop and share key insights.

LABVANTAGE is a thin-client LIMS that is tailored to manage an organization's critical laboratory information across its worldwide R&D pipeline and manufacturing supply chain to optimize productivity and more effectively share knowledge.

REPORTING & ANALYTICS

LABVANTAGE REPORTING & ANALYTICS

Comprehensive Reporting and Analytics-Driven Decision Support

LABVANTAGE offers a number of embedded and integrated reporting and analytics options that can be deployed and securely accessed as needed by the appropriate user. This great variety ensures that LABVANTAGE provides the most versatile and flexible options that deliver reporting and analytics out-of-the-box, while also providing integrative options for reporting and analytics applications you've already invested in.

Combined or deployed as needed, the extensive array of reporting and analytics options offered with LABVANTAGE provide comprehensive, real-time data access, visibility, and reporting capabilities that ultimately arm the organization with the knowledge needed to make more informed critical decisions.

- [LABVANTAGE Dashboards](#) – powerful business management tool provides an “at a glance” view of definable LABVANTAGE pages as well as definable key performance indicators (KPI) using graphical presentation formats including pie charts, spark lines, and gauges.
- [LABVANTAGE Reporting](#) – seamlessly embedded, on-demand reporting and data visualization tools powered by JasperReports.
- [LABVANTAGE Analytics](#) – built-in visual charting for control charts and trend charts powered by JFreeChart.
- [LABVANTAGE Ad-hoc Query Tool](#) – powerful real-time “in-application” data mining and analysis without the need for SQL programming.
- [LABVANTAGE Advanced Reporting](#) – robust information and performance management tools powered by Business Objects.
- [LABVANTAGE Advanced Analytics](#) – full-featured SPC charting and statistics powered by Northwest Analytical.
- [LABVANTAGE Laboratory Intelligence](#) – Business Intelligence for laboratories, multi-source data integration, visualization, and analytics using LIMSVIEW, powered by QlikView.

REPORTING & ANALYTICS


LABVANTAGE REPORTING & ANALYTICS

LABVANTAGE DASHBOARDS

LABVANTAGE Dashboards provide powerful business management tools that offer visual “at a glance” views of definable LABVANTAGE pages as well as definable Key Performance indicators (KPI) using several graphical presentation formats that include trend charts, pie charts, spark lines and gauges. LABVANTAGE Dashboards help you visually ascertain the status and health of your laboratory enterprise by enabling you to:

- Monitor key LABVANTAGE listing pages
- Visually present critical performance measures
- Easily measure efficiencies and inefficiencies across your laboratory enterprise
- Seamlessly monitor your laboratory personnel, equipment and supplies
- Quickly identify negative trends for corrective action
- Highlight new trends online and through detailed reports
- Make more informed decisions based on collected business intelligence
- Drill down directly from the dashboard to the LABVANTAGE data for real time data analysis

LABVANTAGE Dashboards are completely zero-footprint and use no plugs-ins, downloads or applets on the end-user client. Moreover, as an integrated part of LABVANTAGE, it can be tailored to your organization’s end-user requirements via Web Page Designer, the easy-to-use, property-driven configuration tool for LABVANTAGE. Additionally, each user has the ability to personalize their dashboards by picking and choosing which dashboard components to display and use.



Examples of the different dashboard visual indicators such as bar graphs, gauges, trend lines, and pie charts

REPORTING & ANALYTICS

LABVANTAGE REPORTING & ANALYTICS

LABVANTAGE REPORTING

Integrated Reporting Solutions – Powered by JasperReports™

Seamlessly embedded in the application, LABVANTAGE Reports provides on-demand reporting and data visualization. It enables fast design and deployment of both operational and management reports. A number of report templates are included with LABVANTAGE such as Certificate of Analysis, Stability Study Report, Sample Pick Lists, Work Lists, and more.


LABVANTAGE Reports directly supports JasperReports, a 100% Java based solution that uses the JasperReports iReport Designer tool. This reporting application creates a standardized, universal reporting environment, bridging the gap between users in different operational areas. LABVANTAGE significantly simplifies the process of report construction and distribution, helping companies efficiently consolidate and distribute on-demand reports with sample, experiment, and user-specific laboratory data.

Key features & benefits of JasperReports with LABVANTAGE include:

- Leverage's LABVANTAGE's role-based security ensuring proper authentication and access
- Supports LABVANTAGE's multilingual support, including locale formatting, time-zone adjustment, and language translation
- Allows for the embedding of an array of chart types including: tables, cross-tabs, and bar and pie charts, that enhance data visualization
- Ready-to-print report generation delivers prompt data review capabilities in different print and export formats, such as PDF and Microsoft Excel
- Drill-through feature and hypertext enable dynamic linking of reports and supporting documents, along with real-time data analysis
- Ability to version control and deploy multiple report templates, changing the look and feel of the report as desired
- Control and manage amendments and regeneration of reports as results change

Test	Parameter	Result	Unit	Limits	Condition
Physical Properties (Standard)	Appearance	Clear		In: Clear, Transparent, Colorless	Pass
	Odor	None		In: None, Odorless	Pass
	pH	8.4		7.8 - 8.3	Warning

Certificate of Analysis report included with LABVANTAGE



Sample List report generated from LABVANTAGE

REPORTING & ANALYTICS


LABVANTAGE REPORTING & ANALYTICS

LABVANTAGE ANALYTICS

Built-In Control Charts and Trend Charts

LABVANTAGE Analytics provides versatile visual charting. Leveraging LABVANTAGE's integrated trend and control charts, you can effectively and efficiently evaluate processes and initiate process improvements as needed. Control and trend charts provide for the identification and visualization of trends in process data, enabling you to quickly spot problems and highlight critical issues, resulting in the ability to quickly make real time changes to your processes.

Once a chart is displayed, you can interactively zoom in on specific ranges or narrow down their selections to only show a subset of the last number of data points. Graphs can be configured to show violation limits for control, warning, and upper limits to help users identify violations to control limits.



Control and Trend charts displayed in LABVANTAGE

REPORTING & ANALYTICS

LABVANTAGE REPORTING & ANALYTICS

LABVANTAGE AD-HOC QUERY TOOL


LABVANTAGE Ad-hoc Query Tool provides real-time access and visibility to laboratory data. Traditionally, in order to analyze laboratory data, laboratories had to construct multiple sets of reports or queries. Most often these reports or queries are pre-constructed or pre-defined, and managed by a system or database administrator. This can occasionally form a gap between the actual dynamic needs of the end-user and the “canned” information in the report or query. To eliminate this gap and provide users with maximum use and analysis of laboratory data, LABVANTAGE’s Ad-hoc Query Tool enables end-users to construct and execute complex queries against the LABVANTAGE database, with no technical skills, database administration experience, or knowledge of SQL required. Users can execute ad-hoc or saved queries on LABVANTAGE data, including samples, experiments, projects, and results. Built on top of LABVANTAGE’s thin-client enterprise platform, the Ad-hoc Query Tool empowers authorized end-users, both inside and outside of the laboratory, to simply mine the data in LABVANTAGE and extract the information for further use all through a zero-footprint web browser.

End-users can construct their own queries by navigating the object structures in a drag & drop layout, or simply execute previously saved queries by selecting and supplying the criteria needed. Queries can be maintained by the end-user or shared with other users of LABVANTAGE. Moreover, all of the security, compliance, and multilingual support of LABVANTAGE are fully maintained with the Ad-Hoc Query Tool.

LABVANTAGE's Ad-hoc Query Tool empowers end-users to:

- Construct on-demand queries by dragging and dropping appropriate fields onto the query builder palette
- Execute previously saved queries by simply selecting and supplying the criteria required
- Select and move the data columns, as well as sort order on the query results for intuitive viewing
- Export the data to Microsoft Excel for further analysis or processing, as well as, print to PDF
- Access the powerful query interface in any LABVANTAGE list page or through the “Reports” section of LABVANTAGE
- When integrated on a LABVANTAGE list page, users can execute operations on the results such ship, transfer, perform data entry, or add to a work list.

LABVANTAGE’s Ad-hoc Query Tool simplifies access to information and improves knowledge sharing throughout the global enterprise. By expediting flexible queries, organizations gain meaningful results that can be shared to make critical decisions. Ultimately, it gives back the power of query construction to end-users, enabling real-time access and visibility to all data.



An Example of the LABVANTAGE Ad-hoc-Query Interface

REPORTING & ANALYTICS

LABVANTAGE REPORTING & ANALYTICS

LABVANTAGE ADVANCED REPORTING


Advanced Reporting Solutions – Powered by Business Objects™

LABVANTAGE Advanced Reporting provides robust enterprise reporting and analytics. Powered by Business Objects, it offers a complete suite of critical business applications, including information management, performance management, and reporting tools built on the business intelligence platform, Business Objects Enterprise. A robust reporting application, Business Objects offers online reports, as well as, printing, emailing, and exporting of reports to a number of file formats. Moreover, reports can be designed via the thick-client Desktop Intelligence tool, the thin-client Web Intelligence tool, and Crystal Reports.

LABVANTAGE also seamlessly interfaces with Business Objects, enabling enterprises to continue to leverage their existing investment of this enterprise business intelligence tool. Business Objects enables LABVANTAGE users to:

- Leverage the powerful tool set of report authoring, formatting and distribution, enabling rapid creation of flexible, feature-rich reports of LABVANTAGE data.
- Create dashboards and scorecards that summarize and visually represent critical performance and efficiency measures.

LABVANTAGE with Business Objects provides a versatile, interactive viewing framework for companies with a wide range of users, facilitating critical enterprise information sharing.



Examples of Business Objects generated reports.

REPORTING & ANALYTICS

LABVANTAGE REPORTING & ANALYTICS

LABVANTAGE ADVANCED ANALYTICS


SPC Charting and Statistics – Powered by Northwest Analytical® (NWA)

LABVANTAGE Advanced Analytics provides full-featured statistics and charting. Powered by Northwest Analytical, a leading Statistical Process Control (SPC) charting and analysis application, LABVANTAGE offers the standout combination of flexibility and ease of use. It enables a variety of users to graphically analyze and judge the impact of data to make key decisions with minimal training in statistical techniques. Quality Analyst provides comprehensive functionality for companies requiring advanced charting needs for certification, regulatory compliance, process improvement, and cost reduction.

Some of the key benefits that LABVANTAGE offers with NWA Quality Analyst are:

- Multiple chart displays and printing
- Control charts and reports to accompany the Certificates of Analysis (COAs) on any manufacturing process that the laboratory is monitoring
- A state-of-the-art statistical quality/process control (SQC/SPC) charting tool to facilitate the analytical process and discoveries of critical insights through massive loads of data
- Easy access to information needed to foster improved relationships with customers, increase capability, and decrease waste

NWA Quality Analyst is available as an optional charting and analytical application interfaced with LABVANTAGE. Combining a set of comprehensive charting and analysis capabilities, NWA Quality Analyst empowers LABVANTAGE users with a comprehensive choice for their advanced statistics and charting needs.



NWA Control charts displayed in LABVANTAGE

REPORTING & ANALYTICS

LABVANTAGE REPORTING & ANALYTICS

LABVANTAGE LABORATORY INTELLIGENCE

LABVANTAGE offers customers a Business Intelligence (BI) tool called LIMS*View*. Powered by QlikView, LIMS*View* displays market-specific key performance indicators to executives and lab managers, makes ad-hoc and pre-set data analyses that include an endless variety of search criteria, and takes a multisource approach to the BI solution. This means that LIMS*View* can both be integrated with LABVANTAGE and also display data from other enterprise systems. You can combine relevant data into a single BI application, not only the valuable LIMS data, but also information from other systems, like CRM or ERP, data warehouses and spreadsheets.

LIMS*View* helps your business understand its financial performance by examining set targets, profit, volumes and costs. It helps your business understand its throughput by observing each sample status and discovering if the pre-treatment analysis or reporting is your bottleneck. With LIMS*View*, LABVANTAGE offers a solution that relieves the operative system from the burden of continuous business intelligence reporting queries. LIMS*View* extracts and compresses data efficiently to a separate in-memory engine, which serves diverse reporting needs in a dynamic and intuitive way. It also summarizes how your business is doing in production by reviewing quality compliance, trends in analysis and peak results. Ultimately, it helps you identify where you can improve your efforts to be more competitive in your industry. The three major components of LIMS*View* are:

- **QlikView Developer:** Windows-based desktop tool for developing applications
- **QlikView Server:** In-memory analytics engine
- **QlikView Publisher:** Single control point and tool for scheduling, administration and management

Launching LIMS*View* implementation is fast! Your most critical data analysis requirement can be completed in as few as five days. Since LIMS*View* is very scalable, additional requirements can be conveniently solved as needed and can be implemented as a BI tool for deeper analysis and reporting of lab data or create multiple LIMS*View* applications to fulfill various Business Intelligence needs.

SUMMARY

Reporting and analytics plays a vital role in knowledge sharing and informed critical decision-making. LABVANTAGE recognizes that different organizations will have varying requirements for the depth and extent of data accessibility, reporting and analytics. By embedding and integrating multiple leading reporting and analytics applications, LABVANTAGE offers the most versatile and advanced tools required to streamline laboratory management processes and facilitate greater knowledge sharing across your organization and partner network. From the core reporting and analytics functionality, to more intricate, enterprise reporting and analytics capabilities, LABVANTAGE provides the best combination functionality, flexibility, and technology to drive efficiency and knowledge sharing in today's laboratories.



LIMSView powered by QlikView offers a variety of built-in reports and analyses, fully configurable for customer requirements and obtainable as a pre-configured out-of-box solution

LABVANTAGE OVERVIEW

LABVANTAGE Solutions, Inc. is an innovative global provider of enterprise solutions tailored for leading laboratories. We serve numerous laboratories focused on:

- Research & Discovery
- Development (including Life Sciences Pre-clinical & Clinical Development)
- Formulation & Process Research
- Raw Material Testing
- Quality Management

Our global customer base encompasses companies across multiple industries including:

- Pharmaceutical & Biosciences
- Health & Personal Care
- Food & Beverage
- Process Chemicals
- Petrochemicals
- Parts, Hardware & Equipment
- Law Enforcement & National Security
- Contract Services

LABVANTAGE is a thin-client Laboratory Information Management Solution (LIMS) tailored to manage an organization's critical laboratory information across its worldwide R&D pipeline and manufacturing supply chain in order to optimize productivity and more effectively share knowledge.

Headquartered in Somerset, New Jersey, United States, LABVANTAGE also has offices in High Wycombe, United Kingdom, Monza, Italy, and Kolkata, India, and utilizes a strategic network of global business partners to provide regional sales and customer support.

LABVANTAGE is a privately held affiliate of The Chatterjee Management Group Company, a \$1.5 billion private investment group that invests globally in infrastructure, communications and information technology companies.

LABVANTAGE Brochure

WWW.LABVANTAGE.COM

©2013 LABVANTAGE Solutions, Inc. All rights reserved.

For information regarding this and other LABVANTAGE products, please contact LABVANTAGE Solutions, Inc. The information contained in this document is proprietary and confidential to LABVANTAGE Solutions, Inc. No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording, for any purpose, without the express written permission of LABVANTAGE Solutions, Inc.

This document is subject to change without notice, and LABVANTAGE does not warrant that the material contained in this document is error-free. If you find any problems with this document, please report them to LABVANTAGE in writing. LABVANTAGE® and Web Page Designer are trademarks of LABVANTAGE Solutions, Inc. Red Hat® and JBoss® are registered trademarks of Red Hat, Inc. Microsoft®, SQL Server and Internet Explorer are either registered trademarks or trademarks of Microsoft Corporation. Apple Safari® is a registered trademark of Apple Inc. Oracle®, WebLogic, and Java are registered trademarks of Oracle and/or its affiliates. IBM® and WebSphere are trademarks of International Business Machines Corporation. Jaspersoft® and JasperReports are trademark of Jaspersoft Corporation. SAP® BusinessObjects is the trademark or registered trademark(s) of SAP AG. Intel® and Xeon are trademarks of Intel Corporation. Sybase Enterprise Application Server is a Obsolete Sybase Mark. All other trademarks are property of their respective owners.

This document may contain statements concerning possible functionality for LABVANTAGE software products and technology. LABVANTAGE disclaims any express or implied commitment to deliver functionality or software unless or until generally commercially available. Any statements of possible future direction are for information purposes only, and LABVANTAGE makes no express or implied commitments or representations concerning the timing and content of any future functionality or releases.